

Ingenuity in
infection control
since 1830



AIRGENE CE
disperses
disinfectant to
difficult &
inaccessible
surfaces.
Open your
practice with
confidence.

AIRGENE CE

**Exclusive
to Eschmann**

Protecting your patients, your staff and your business from the risk of infection from harmful pathogens is of paramount importance, especially in today's world.

AIRGENE CE is for professional use and designed for rapid decontamination of inaccessible surfaces to reduce cross infection in high risk environments.

AIRGENE CE

Surgery Disinfection

AIRGENE CE Airborne Disinfectant has been designed to disperse high-level disinfectant, decontaminating surfaces and hard to reach areas, reducing the risk of cross infection in high risk environments.

Once the valve is released, the AIRGENE CE forms a cloud in approximately 1 minute, covering up to 40m³.

AIRGENE CE takes just one hour to achieve a double effect: decontamination of micro-organisms on all surfaces and hard to reach areas and the effective elimination of unpleasant odours caused by bacterial decomposition.

The AIRGENE CE is to be used in conjunction with your existing cleaning protocols.

- Unique design allows for dispersion of high-level disinfectant to difficult and inaccessible places.
- In areas where the risk of micro-particle contamination is especially high (e.g high-speed handpieces in a dental surgery can spread particles over a wide area), the daily use of AIRGENE CE is recommended.
- 50ml canister covers an area of up to 40m³ (4 x 4 x 2.5m).
- Place on a flat surface in the centre of a room - release the valve - leave for just one hour.
- Free of flammable and toxic ingredients.

AIRGENE CE is registered as a class IIa medical device, CE marked for use in healthcare. It is widely effective against bacteria, fungi and viruses as demonstrated in independent tests carried out under EN protocols:

BACTERICIDAL:

- EN 1276 - Gram+ (Staphylococcus aureus, Enterococcus hirae)
Gram- (Pseudomonas aeruginosa, Escherichia coli)
- EN 13697** - Gram+ (Staphylococcus aureus, Enterococcus hirae).
Gram- (Pseudomonas aeruginosa, Escherichia coli)
- EN 13727 - Gram+ (Staphylococcus aureus, Enterococcus hirae)
Gram- (Pseudomonas aeruginosa)

MYCOBACTERICIDAL:

- EN 14348 - Mycobacterium avium, Mycobacterium terrae

FUNGICIDAL:

- EN1650 - Candida albicans, Aspergillus niger
- EN 13697** - Candida albicans, Aspergillus niger
- EN 13624 - Candida albicans, Aspergillus brasiliensis

VIRUCIDAL:

- EN 14476 - Influenza A (H1N1)*

*Surrogated flu virus for lipophilic viruses (Flu, Coronavirus, Ebola, Hepatitis, HIV...)

**Product tested without alcohol

Pack options and prices

Cat No.	Part Description	Price (excl.VAT)
AG012CE	AIRGENE CE (single unit)	£14.20
AG012CEX12	AIRGENE CE (Case of 12)	£156.20

Coverage



Eschmann Technologies Limited
Eschmann House, 15 Peter Road, Lancing,
West Sussex, BN15 8TJ United Kingdom

t: +44 (0)1903 753322
e: info@eschmann.co.uk
w: eschmann.co.uk

Ingenuity in
infection control
since 1830

