
Ingenuity in infection control since 1830





# Little Sister UltraSONICA

Exclusive to Eschmann

The Little Sister Ultrasonic Cleaners utilise the advanced Sweep System technology to provide effective and reliable cleaning of instruments in the shortest possible time.

With a Little Sister UltraSONICA, you can receive annual testing with validation and User Training, combined with unrivalled, nationwide support from Eschmann.

## Little Sister UltraSONICA

### The Little Sister UltraSONICA uses the latest advanced Sweep System technology for faster and improved cleaning of instruments.

This ultrasonic technology uses a process known as cavitation to safely and effectively clean instruments. This process combined with the Hands Free Instrument Transfer (HFiT) Plus Cassettes, validated for use with the UltraSONICA range, also reduces the risk of injury and cross infection which may occur during manual cleaning.

Available in 3 sizes to suit different processing needs, the Little Sister UltraSONICA incorporates the latest oscillator and transducer technologies with simple, easy to use controls, offering faster and better cleaning than more traditional ultrasonic cleaners.

#### **Consumables & Accessories**



Cleaning Test Strips x200 Part number: 010.008.0003

Support Base for Cleaning Test x1 Product number: 090.004.0072

Essential Protein Test Pen x25
Product number: 3371



**Log Book for Ultrasonic Cleaners** Product number: 119878

Stainless Steel Container for Small Parts Product number: 090.004.0074



**HFiT Plus Instrument Cassettes** 

Note: Silicones available in 6 colours. Product code last digit 'x' determines the colour. Replace 'x' with: 2-Blue, 3-Red, 4-Green, 5-Yellow, 6 Silver or 8-White. Product number: 872280x (HFiT-5), 872281x (HFiT-10), 872282x (HFiT-20)

### **Specifications & Prices**

Little Sister UltraSONICA	2200m Ref: 090.003.0002	2400m Ref: 090.003.0006	3200m Ref: 090.003.0010
Tank Capacity (litres)	3.0	4.5	6.0
Transducers	2	2	3
Useable Chamber Size mm (w x l x d)	240 x 140 x 100	300 x 150 x 100	300 x 240 x 100
Overall Dimensions mm (w x l x h)	270 x 170 x 260	325 x 175 x 260	325 x 270 x 260
Weight kg (empty)	3.2	3.5	4.3
Accessories Included (x1 of each)	Stainless Steel Cover Stainless Steel Basket	Stainless Steel Cover Stainless Steel Basket	Stainless Steel Cover Stainless Steel Basket

#### **Eschmann Technologies Limited**

Eschmann House, 15 Peter Road, Lancing, West Sussex, BN15 8TJ United Kingdom

t: +44 (0)1903 753322 e: info@eschmann.co.uk w: eschmann.co.uk Ingenuity in infection control since 1830

