

Miele

High quality **laundry and washer disinfector solutions** for veterinary practices

from Eschmann and Miele

Quality is at the heart of everything Eschmann and Miele provide. All Miele products are designed, built and tested to operate in a busy commercial environment for a minimum of ten years. For example, for our laundry washers that's 30,000 cycles, which equates to 8 cycles a day, every day, for ten years. You would typically replace a domestic washing machine 6 times in that period for the same use.

Combined with Care & Cover service offering from Eschmann and Miele, you can have peace of mind that your equipment is compliant, running with optimal efficiency and minimal downtime. Care & Cover provides all the technical support you need through our nationwide team of engineers, and includes unlimited breakdown cover, call-outs, parts and labour.



Washer disinfector features:

- ✓ Reliability and longevity
- ✓ Fast wash cycle times
- ✓ WRAS approved
- ✓ Large capacity, can fit up to 360 instruments
- ✓ Integrated water softener
- ✓ Quiet operation

Laundry washers and dryers features:

- ✓ Reliability and longevity
- ✓ Tested to 30,000 wash cycles
- ✓ Fast cycle times
- ✓ WRAS approved
- ✓ Stackable
- ✓ Easy to use

Specifications: PG8581 washer disinfector



PG8581
Washer disinfector

With each PG8581 washer disinfector for veterinary practices, we include:

- 1 x A202 lower basket
- 1 x A103 upper basket
- 2 x E142 half width mesh tray


PLUS

- 1 x detergent (5L), 1 x rinse aid (5L), 1 x pellet salt (10kg), 2 x additional instrument baskets and 1 x washer disinfector log book.

Cycle times (mins.)	From 50 mins*
Number of programmes	3
Water connections	Cold water
External dimensions H/W/D (mm)	835/600/600
Steam condenser	•
Integrated water softener	•
Electrical connection 1N AC 230V 50Hz	•
Total rated load (kW)	6.3
Electrical connection 3N AC 400V 50Hz	•
Total rated load (kW)	9.3

* Without drying, depending on the load


Accessories for the PG8581



A103 upper basket



E142 half width mesh tray



A202 lower basket

Specifications: Little Giant laundry washers and dryers



Evolution PWM 307 [EL DP]



Little Giant PWM 507 [EL DV]



Evolution PDR 307 HP [EL]



Evolution PDR 307 [EL]

WRAS	Cat 5	Cat 5
Honeycomb drum with perforated drum (rear panel)	•	•
Programme Packages	Household	Full suite ¹
Programme duration in minutes with cold water connection	59	59
Programme duration in minutes with cold and hot water connection	-	49
Hot water connection	-	•
Spin Speed	1,600	1,600
g factor	704	704
Tested operating hours	30,000	30,000
No. of cycles per day	Unlimited	Unlimited
Drainage type	Drain pump	Drain valve

Drum volume in l	130	130
Maximum evaporation capacity in l/hour	2.7	2.6
Energy consumption for Normal programme in kWh/kg	0.17	0.49
Programme Packages	Household	Household
Programme duration on Normal mixed fabrics (40% RM) in min	57	66
Programmability	Programmable	Programmable
Tested operating hours	30,000	30,000
Energy efficiency class	A++ (higher power motor)	C
Drying system	Heat pump	Vented

¹ - includes hygiene cycle which is compliant to HTM0104

Why Eschmann and Miele

Eschmann and Miele have worked in partnership for 20 years. The expertise and unrivalled service offering from Eschmann combines perfectly with the high quality and longevity of Miele products, providing reliable high quality solutions to our customers. To ensure complete peace of mind, we also offer our customers a range of preventative service and maintenance contracts together with a first class emergency repair response.

Eschmann Technologies Limited
Eschmann House, 15 Peter Road, Lancing, West Sussex, BN15 8TJ United Kingdom
t: +44 (0)1903 753322 e: info@eschmann.co.uk w: eschmann.co.uk

